



www.pcultrasound.com
info@pcultrasound.com



TELEMED
ULTRASOUND MEDICAL SYSTEMS



WHY TELEMED'S ARTUS COMPACT?

- Powered via USB C/USB3.0
- Small footprint and high-performance
- A broad range of transducers for various applications
- Attractive price

AFFORDABLE YET POWERFUL

ArtUs C is a portable USB-powered ultrasound system developed to fill the gap between ArtUs and MicrUs scanners. Offering the same technologies as its larger brother, ArtUs C(ompact) is priced at half as much as its ArtUs. It is powered by a USB-C port, eliminating the need for external power sources, and equipped with Color, Power, Directional Power, Pulsed Wave Doppler, Spatial Compound, Harmonic Imaging. Optional RF Data transfer and Extensive triggering capabilities make it an attractive tool for universities and research labs.

Convex array transducer

C5-2N60-D3 Radius 60mm 2-5MHz
For abdominal, obstetrics, urology, veterinary

Linear array transducers

L12-5N25-D3 Length 25mm
L12-5N40-D3 Length 40mm
L12-5N50-D3 Length 50mm 5-12MHz
For MSK, vascular, small parts, veterinary

Phased array (Sector) transducer

P4-2R20-D3 2-4MHz
For cardiology, adult

Microconvex and Endocavity transducers
coming soon

CALL US: +370 5 2106272