



www.pcultrasound.com
info@pcultrasound.com



TELEMED
ULTRASOUND MEDICAL SYSTEMS



WHY TELEMED'S SPECTRUS?

- Three-In-One: linear, convex and sector scanning
- One of the most light-weight probes on the market: 220g
- Optimised power consumption: 3 hours scanning time
- Wireless smart charger/case
Windows, Android and iOS compatible

TWO HEADS MULTIPLE APPLICATIONS

SpectrUs features a compact, lightweight, and portable design with extended battery life and wireless connectivity. Its operating modes include B-mode, Compound Imaging, Harmonic Imaging, SRI (Speckle Reduction Imaging), Color Flow, Power Doppler, Directional Power Doppler, and an optional Needle Enhance mode.

Equipped with a sensitive, broadband convex transducer, the SpectrUs is capable of performing cardiac sector imaging, supporting cardiology applications.

Convex array transducer

For abdominal, obstetrics, urology, paediatrics, cardiology, veterinary

Linear array transducer

For vascular, MSK, lungs, paediatrics, small parts, veterinary

CALL US: +370 5 2106272