

MICRUS PRO



TELEMED




ULTRASOUND SCANNERS

Located in Vilnius, Lithuania, TELEMED designs, develops and manufactures PC-based ultrasound scanners and related equipment. Using modern hardware and software technologies, TELEMED offers cost-effective and powerful solutions to clinics, private doctors and original equipment manufacturers.

TELEMED was established in 1992 for research and design in the field of ultrasound medical equipment. Since 1995, when first PC-based scanner was shipped to the customer, easy-to-use and scalable TELEMED's solutions have become an important research tool not only for physicians, but also for scientific institutions.

Version: A

 OFFICE ADDRESS: Savanoriu pr. 178A, VILNIUS, LT-03154, LITHUANIA

 PHONE NUMBER: (+370-5) 2106272

 EMAIL ADDRESS: info@telemed.lt

 WEBSITE: www.pcultrasound.com

MICRUS PRO IS A HANDHELD PERSONAL ULTRASOUND SCANNER BASED ON A POPULAR MICRUS PLATFORM, ITS UPGRADED ELECTRONICS AND CAREFULLY DESIGNED PRINTED CIRCUIT BOARDS PROVIDE IMAGE QUALITY APPRECIATED BY MICRUS CUSTOMERS WHILE THE SIZE IS MINIMIZED TO FIT INTO A TRANSDUCER HOUSING.

MICRUS PRO IS AN APPLICATION-DRIVEN DEVICE AND WORKS WITH WINDOWS AND ANDROID PLATFORMS ON MODERN TABLETS AND MOBILE PHONES.

SOFTWARE DEVELOPMENT KIT (SDK) FOR WINDOWS AND ANDROID CAN BE PROVIDED ON REQUEST FOR INDUSTRY, DEVELOPERS OR RESEARCH TEAMS.



MICRUS PRO-L40S

- B/W ULTRASOUND SCANNER INTENDED FOR SCREENING PURPOSES
- POWERED VIA USB INTERFACE
- 2.5 HOURS CONTINUOUS WORK
- REMOVABLE USB CABLE AND CAP FOR EASY DISINFECTION
- MS WINDOWS AND ANDROID SOFTWARE AVAILABLE
- 5-12MHZ 40MM ULTRASOUND TRANSDUCER



MICRUS PRO-L40N

- B/W ULTRASOUND SCANNER INTENDED FOR SCREENING PURPOSES
- POWERED VIA USB INTERFACE
- 2.5 HOURS CONTINUOUS WORK
- REMOVABLE USB CABLE AND CAP FOR EASY DISINFECTION
- MS WINDOWS AND ANDROID SOFTWARE AVAILABLE
- 5-12MHZ 40MM ULTRASOUND TRANSDUCER



MICRUS PRO-C60S

- B/W ULTRASOUND SCANNER INTENDED FOR SCREENING PURPOSES
- POWERED VIA USB INTERFACE
- 2.5 HOURS CONTINUOUS WORK
- REMOVABLE USB CABLE AND CAP FOR EASY DISINFECTION
- MS WINDOWS AND ANDROID SOFTWARE AVAILABLE
- 2-5MHZ R60MM ULTRASOUND TRANSDUCER