



Compatible Disinfectants by TELEMED Transducer Model

Rev.B, Effective date 2021-05-24

TELEMED Transducer	Cleaners	Low Level Disinfectants	High Level Disinfectants
BIPC6.5/10/128Z-4	Klenzyme		Cidex OPA, Nu-Cidex, Salvianis pH7/pH10
BIPC8-4R10N-4	Klenzyme		Cidex, Cidex Plus, Cidex OPA, Nu-Cidex, Salvianis pH7/pH10
BIPC9-4R10H-4	Klenzyme		Cidex OPA, Nu-Cidex, Salvianis pH7/pH10
BIPL10-4L60N-4	Klenzyme		Cidex, Cidex Plus, Cidex OPA, Nu-Cidex, Salvianis pH7/pH10
BIPL12-5L70H-4	Klenzyme		Cidex OPA, Nu-Cidex, Salvianis pH7/pH10
BIPL7.0/60/128Z-4	Klenzyme		Cidex OPA, Nu-Cidex, Salvianis pH7/pH10
C3.5/20/128Z-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
C3.5/20/64D-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
C3.5/60/128Z-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
C3.5/60/64D-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
C4.5/50/128Z-2	Cidezyme / Enzol		Cidex OPA
C5-2H60-A5	Cidezyme / Enzol, MetriZyme	Milton antibacterial solution	Cidex, Cidex Plus, Cidex OPA
C5-2R60HI-5	Cidezyme / Enzol, MetriZyme	Milton antibacterial solution	Cidex, Cidex Plus, Cidex OPA
C5-2R60NI-5	Cidezyme / Enzol, MetriZyme	Milton antibacterial solution	Cidex, Cidex Plus, Cidex OPA
C5-2R60S-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
C6-1H50-A5	Cidezyme / Enzol, MetriZyme	Milton antibacterial solution	Cidex, Cidex Plus, Cidex OPA
C7-3R50NI-5	Cidezyme / Enzol, MetriZyme	Milton antibacterial solution	Cidex, Cidex Plus, Cidex OPA
EC6.5/10/128Z-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
EC6.5/10/64D-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
HL9.0/40/128Z-4	Klenzyme		Cidex, Cidex Plus, Cidex OPA, Nu-Cidex, Salvianis pH7/pH10
HL9.0/40/64D-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
HL9.0/60/128Z-2	Cidezyme / Enzol		Cidex OPA
L12-5L40N-4	Klenzyme		Cidex, Cidex Plus, Cidex OPA, Nu-Cidex, Salvianis pH7/pH10
L12-5L40S-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
L12-5N40-M3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
L15-6L25N-4	Klenzyme		Cidex, Cidex Plus, Cidex OPA, Nu-Cidex, Salvianis pH7/pH10
L15-6L25S-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
L15-7H40-A5	Cidezyme / Enzol, MetriZyme	Milton antibacterial solution	Cidex, Cidex Plus, Cidex OPA
L15-7L40H-5	Cidezyme / Enzol, MetriZyme	Milton antibacterial solution	Cidex, Cidex Plus, Cidex OPA
L18-10L30H-4	Klenzyme		Cidex, Cidex Plus, Cidex OPA, Nu-Cidex, Salvianis pH7/pH10
L18-7H30-A5	Cidezyme / Enzol, MetriZyme	Milton antibacterial solution	Cidex, Cidex Plus, Cidex OPA
LF11-5H60-A3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
LF9-5N60-A3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
LV7.5/60/128Z-2	Cidezyme / Enzol		Cidex OPA
LV7.5/65/64D-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
LV8-4L65S-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
LV8-5L60N-2	Cidezyme / Enzol		Cidex OPA
LV8-5N60-A2	Cidezyme / Enzol		Cidex OPA
MC10-4R12N-4	Klenzyme		Cidex, Cidex Plus, Cidex OPA, Nu-Cidex, Salvianis pH7/pH10
MC10-5R10S-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
MC4-2R20S-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
MC8-4R20S-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
MCV9-5N10-A3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
MCV9-5R10N-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
MCV9-5R10S-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
P5-1L15SI-6	Klenzyme, Transeptic spray, MetriZyme, Cidezyme / Enzol	Super Sani-cloth, T-Spray	Anioxyde 1000, Cidex, Cidex Plus, Cidex OPA, Gigasept FF, Sporox II, Virkon S, Wavicide-01
P5-1S15-A6	Klenzyme, Transeptic spray, MetriZyme, Cidezyme / Enzol	Super Sani-cloth, T-Spray	Anioxyde 1000, Cidex, Cidex Plus, Cidex OPA, Gigasept FF, Sporox II, Virkon S, Wavicide-01
P8-3L10SI-6	Klenzyme, Transeptic spray, MetriZyme, Cidezyme / Enzol	Super Sani-cloth, T-Spray	Anioxyde 1000, Cidex, Cidex Plus, Cidex OPA, Gigasept FF, Sporox II, Virkon S, Wavicide-01
PV6.5/10/128Z-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA
PV6.5/10/64D-3	Incidin OxyFoam	Incidin OxyFoam	Cidex, Cidex Plus, Cidex OPA